

Table Top Autoclave

APPLICATIONS FOR

Aureole Table Top High Pressure High Vacuum High Speed Sterilizers are ideal for all applications requiring rapid, reliable and total destruction of all types of living Microorganisms. Ideal for Research lab, Hospitals, Dental clinics, Operation theaters, QC/Micro etc.



Table Top Autoclave

Model Available are in liters: JE-T- 8 / 12 / 16 / 18 / 23. (B- CLASS & N CLASS)

FEATURES & ADVANTAGES

- Automatic computer control, LCD dynamic display, fault self - diagnosis, display the error code.
- Reference European class B standard, three pre vacuum and vacuum deep drying function, to make the sterilization and drying sterilization, drying to achieve the best working conditions.
- Preset sterilization procedures, sterilization and the drying time can be adjusted.
- SS 316 L stainless steel containers and a sealing cover, double pressure protection safety door lock.
- Start working procedures prevent when door is not locked well.
- Safety pressure protection lock, automatic induction sterilizing cavity pressure, can't open door during sterilization process.
- Built -in independent steam generator, fast saturated steam, so that the high temperature steam fully reaches every corner.
- Simple and convenient operating interfaces, LCD dynamic display, display the status information in real time.
- Open type top water tank, filling and cleaning for user.
- With liquid sensor for water level, to remind user to add water or drainage.
- The built-in safety valve, automatically exhaust and unloading pressure when over pressure.
- Equipped with temperature protection device, in the cavity wall temperature and steam generator device is too high will automatically disconnect.
- With B&D test, Helix Test, and Vacuum test program.
- With optional external printer.
- The optional USB function and U disk, the data stored in the U disk, connect the computer to view and save.



Aureole Process Equipment Pvt. Ltd.

www.aureoleproequip.com